

# D-DIVIDERMAKS

SAMPLE DIVIDER



D-Dividermaks model sample divider (Riffle) is used to mix product samples such as wheat, corn, rye, oats, paddy, rice, lentils etc. quickly, homogeneously and standardly.

- Ability to mix samples homogeneously
- Ergonomic design suitable for desktop use
- Provides fast and reliable results
- 18 channels available for homogeneous mixing
- Collection containers have a capacity of 1 kg

**MAKSWELL**  
innovation

MAKSWELL MAKİNE GIDA LTD. ŞTİ  
İvedik O.S.B. Mah. 1514 Cad. No:3  
Yen mahalle / ANKARA / TURKEY

Tel: +90 (312) 395 67 37  
Mail: [sales@makswell.com.tr](mailto:sales@makswell.com.tr)  
Web: [www.makswell.com.tr](http://www.makswell.com.tr)