

**Digital Control Unit**

Temperature can be programmed between 0-100 °C with 3-digit LED display control unit

**Integrated Work**

When the water temperature in the boiler of the Dual Enzymaks device exceeds 92 degrees, Coolmaks automatically starts to work

**Double Fan**

Equipped with two fans with automatic speed calibration for efficient cooling

**Ergonomic Design**

With its ergonomic design, it saves space in the laboratory and can be mounted behind the Dual enzymaks.

During the falling number analysis at boiling temperature, the bath lid must be constantly cooled to prevent water loss. The Coolmaks device cools the bath lid by circulating 500 ml of water in the pump tank. The device works with a circulation system, preventing water waste and providing efficient cooling.

**MAKSWELL**  
innovation

MAKSWELL MAKİNE GIDA LTD. ŞTİ  
İvedik O.S.B. Mah. 1514 Cad. No:3  
Yenimahalle / ANKARA / TURKEY

Tel: +90 (312) 395 67 37  
Mail: sales@makswell.com.tr  
Web: www.makswell.com.tr