



Digital Control Unit

The time can be set between 0 and 99 minutes with the 4-digit LED display control unit

Teflon Coated Plates

The surfaces of the plates are coated with Teflon to minimize surface tension

Double Thermostat

Equipped with two thermostats for 100°C standby and 150°C test temperature

Latch Lock

Used to keep the distance between two plates constant

Dry gluten device (Drymaks) is used to dry wet gluten. Wet gluten samples placed between two plates at 150 ° C are dried for 4 minutes to calculate the water absorption and dry gluten value of the gluten. Drymaks has a two-sample capacity and saves on personnel costs and time.

MAKSWELL
innovation

MAKSWELL MAKİNE GIDA LTD. ŞTİ
İvedik O.S.B. Mah. 1514 Cad. No:3
Yen mahalle / ANKARA / TURKEY

Tel: +90 (312) 395 67 37
Mail: sales@makswell.com.tr
Web: www.makswell.com.tr