

Accessories supplied with Sedimaks device

List of required accessories for sedimentation analysis is as below. All the accessories may be obtained from our company after the sale.



No	Description	Quantity
----	-------------	----------

1	Reservoir Bottle, 1 Lt	2
2	Automatic Pipette, 25 ml	1
3	Automatic Pipette, 50 ml	1
4	Lactic Acid Solution, 1 Lt	1
5	Bromine Pheno Solution, 1 Lt	2
6	Anlysis Tube	6
7	Plastic Funnel	1
8	Electric Power Cable	1
9	Touch screen Pen	1
10	Cleaning Brush	1

Recommended Equipment

Precision Scale (0.01 sensitive)



Accurate measurement of sample quantity is an important factor in terms of analysis results. Therefore, reliable and appropriately sensitive scales should be used. Our company provides scales to our customers upon request.